













## SUMMARY OF TRANSPORT MEDIA FOR MICROBIOLOGY


Please refer to Laboratory Service Guide for further details.







<u>Specimen type</u>	<u>Test required</u>	<u>Container</u>
Blood & Bone Marrow	<b>Aerobic culture</b>	 <b>Bact/Alert bottles (▲)</b> <u>Aerobic</u> - Adult (green - FA) - Paeds (yellow - PF) <u>Anaerobic</u> - Orange (FN)
	<b>Aerobic &amp; Anaerobic culture</b>	
	<b>Fungus culture (Candidasis only)</b>	
	<b>Fungus culture (for dimorphic fungi and mould)</b>	 Heparin tube (▲)
	<b>Acid Fast Bacilli (AFB) culture</b>	 Heparin tube (▲) (EDTA tube not suitable)
	<b>Molecular Detection of <i>Mycobacterium tuberculosis</i></b>	 EDTA tube (▲)
Blood (serum)	<b>ASOT</b>	 Plain (red top) or Gold top tube (▲)
	<b>Mycoplasma Total Antibody</b>	
	<b>Aspergillus galactomannan Antigen</b>	
	<b>Cryptococcal Antigen</b>	
Body Fluid, Body Aspirate	<b>Aerobic &amp; Anaerobic culture</b>	 Sterile container (▲)
	<b>Fungus culture</b>	
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>	






CSF	<b>Aerobic culture</b>	 <p>Sterile container (▲)</p>
	<b>Fungus culture</b>	
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>	
	<b>Cryptococcal Antigen (CSF)</b>	
	<b>Aspergillus galactomannan Antigen (CSF)</b>	
	<b>Meningitis and Encephalitis Panel PCR</b>	



Eye / Cornea Scraping	<b>Aerobic culture</b>	 <p>Chocolate plate; CDC ANA Blood plate (ㄨ)</p>
	<b>Aerobic &amp; Anaerobic culture</b>	
	<b>Fungus culture</b>	 <p>Sabouraud's plate with chloramphenicol (ㄨ)</p>
	<b>Acanthamoeba culture</b>	 <p>Page Saline (ㄨ)</p>
	<b>Microsporidia microscopy (eye / cornea)</b>	 <p>Slide (ㄨ)</p>


Genital tract swab	<b>Aerobic culture</b>	 <p>eSwab (Single swab) (▲)</p>
	<b>Trichomonas / Candida / Bacterial vaginosis</b>	
	<b>Streptococcus Group B screen</b>	



Wound and miscellaneous site swab	<b>Aerobic culture</b>	 <p>eSwab (Single swab) (▲)</p>
	<b>Aerobic &amp; Anaerobic culture</b>	
	<b>Fungus culture</b>	
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>	


Respiratory specimens (sputum, ETT aspirate, BAL), nasal and throat swabs	<b>Aerobic culture</b>	 Sterile container (▲)   eSwab (Single swab) (▲) (for Nasal & Throat swab)
	<b>Fungus culture</b>	
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>	
	<b>Molecular Detection of <i>Mycobacterium tuberculosis</i></b>	
	<b>Aspergillus galactomannan Antigen (BAL)</b>	
Respiratory: Nasopharyngeal swab / Bronchoalveolar lavage (BAL)	<b>Respiratory Pathogens Panel PCR</b>	 Universal transport medium (UTM) and Flocked swab (ꠞ)   Sterile container (▲)
Respiratory: Nasopharyngeal swab	<b>Influenza A, B and RSV PCR</b>	 Universal transport medium (UTM) and Flocked swab (ꠞ)
Skin scraping, nail clipping, hair	<b>Fungus smear</b>	 Sterile container (▲)
	<b>Fungus culture</b>	
	<b>Scabies microscopy (skin scraping)</b>	

Nasal, Axilla, Groin swabs	<b>MRSA screen</b>	 MRSA eSwab (Double swabs) (▲)
	<b><i>Candida auris</i> DNA PCR</b>	
	<b><i>Candida auris</i> screen</b>	
Rectal swab	<b>CRE screen</b>	 eSwab (Single swab) (▲)
	<b>VRE screen</b>	
Stool	<b>Aerobic culture</b>	 Sterile container with spatula (▲)
	<b>Occult blood</b>	
	<b>Ova, cyst and parasites examination</b>	
	<b><i>C. difficile</i> toxin</b>	
	<b>Faecal Fat</b>	
	<b>Cryptosporidium, Cyclospora, Cystoisospora Microscopy</b>	
	<b>Gastrointestinal Pathogens Panel PCR</b>	
	<b>Microsporidia Microscopy (stool)</b>	
	<b>CRE screen (stool)</b>	
	<b>VRE screen (stool)</b>	
	<b>Faecal Calprotectin</b>	
Tissue / Biopsy / Bone	<b>Aerobic &amp; Anaerobic culture</b>	 Sterile container (▲)  Gastric biopsy for <i>H. pylori</i> culture: send biopsy in 2 – 3 drops sterile saline. The container should be kept upright, tissue must be submerged at all times.
	<b>Fungus culture</b>	
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>	
	<b>Molecular Detection of <i>Mycobacterium tuberculosis</i></b>	
	<b><i>Helicobacter pylori</i> culture (gastric biopsy)</b>	
Urethral, Cervical, Neonatal eye	<b>Gonococcus culture</b>	 eSwab (Single swab) (▲)

Urine			Boric acid container (Adult) (▲)
	<b>Aerobic culture</b>		Boric acid container (Paeds) (ㄸ)

Urine	<b>Fungus culture</b>		Sterile container (▲)
	<b>Acid Fast Bacilli (AFB) culture (Volume ≥ 40mL)</b>		
	<b>Molecular Detection of <i>Mycobacterium tuberculosis</i></b>		
	<b>Urine Red Cell Morphology</b>		

Other Microbiology specimens (please specify collection site)	<b>Aerobic culture</b>		Sterile container (▲)
	<b>Aerobic &amp; anaerobic culture</b>		
	<b>Fungus culture</b>		
	<b>Acid Fast Bacilli (AFB) culture &amp; smear</b>		
	<b>Molecular Detection of <i>Mycobacterium tuberculosis</i></b>		
	<b>Gram smear</b>		
	<b>Fungus smear</b>		
	<b>AFB smear</b>		
			eSwab (Single swab) (▲)

Various cultures and smears	<b>Amies Transystem gel swabs are still applicable for non eSwab users from non-NUH institutions</b>		Amies Transystem (*)
-----------------------------	--	--	----------------------

(▲): For in-house NUH patients, please obtain from MMD ext 25196.

(ㄸ): Please obtain container from NUH Microbiology Laboratory ext 24343.

Note: For non-NUH clients, please liaise with National University Hospital Referral Laboratory (65) 6778 5171.